Digital Thermometer English

A. Sensor Tip
B. Display Window
C. I/O Button
D. Storage Case
Introduction
Thank you for purchasing the OMRON Digital Thermometer Eco Temp Basic.

Important Safety Notes
To assure the correct use of the product basic safety measures should always be followed including the precautions listed below.
Warning
A high or prolonged fever requires medical attention, especially for small children.

Table with 3 columns: Indicator, Cause, and Remedy. Rows include Low (L), High (H), and Air Temperature Displayed.

7. Storage
Store the thermometer in its storage case.
Do not store the thermometer in the following types of places. Doing so may damage the thermometer.

8. Cleaning and Disinfecting
OMRON recommends to wipe the device with damp cloth. For disinfection, 70% Ethanol or Isopropyl alcohol can be used.

9. Specifications
Product Type and Number: Digital Thermometer MC-246-EMC-246-EB
Product Name: Eco Temp Basic
Power Supply: 1.5V DC, 1LR41 Alkaline-Manganese Button Battery

Applied part: Type B
Notes:
The specifications may be changed without prior notice.
This unit is water-resistant.

Digital thermometer Model Eco Temp Basic

- Digitales Fieberthermometer
Thermomètre digital
Termometro digitale
Termómetro digital
Digitalt termometer
Digitaalinen kuumemittari
Digitalt feberthermometeret
Digitale koortsthermometer
Ψηφιακό θερμομέτρο
Dijital klinik termometresi
Цифровой термометр

Large multilingual table with columns for various languages including English, Spanish, French, Italian, Portuguese, Danish, Finnish, Norwegian, Swedish, Dutch, Greek, Turkish, Russian, and Arabic. Each column contains instructions for using the thermometer in that language.

10. Display Indications and Troubleshooting
This product is guaranteed by OMRON for 3 years after date of purchase. Any damage caused by improper handling shall not be covered by warranty.

Read the Instructions in this manual carefully.
CE0197

11. Battery Replacement
Battery: Alkaline-Manganese Button Battery LR41 (available commercially)
Battery Replacement Indicator (BA)
Replace the battery when the battery replacement indicator appears when the thermometer is switched on.

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OMRON HEALTHCARE UK LTD. Opal Drive, Fox Milne, Milton Keynes, MK15 0DG, U.K.